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
1. Sissener T. Suture Patterns. Comp Anim. 2006;11:14-19. 2. Wong N. Review of Continuous Sutures in Dermatologic Surgery. J Dermatol Surg Oncol. 1993;19:923-931. 3. Boutros S, Weinfeld A, Friedman J. Continuous Versus Interrupted Suturing of Traumatic Lacerations: A Time, Cost, and Complication Rate Comparison J Trauma Injury Infect Crit Care. 2000;48:495-497. 4. Seiler C, Bruckner T, Diener M et al. Interrupted or Continuous Slowly Absorbable Sutures For Closure of Primary Elective Midline Abdominal Incisions A Multicenter Randomized Trial. Ann Surg. 2009;249:576-582. 5. Gupta H, Srivastava A, Menon G, Agrawal C, Chumber S, Kumar S. Comparison of Interrupted Versus Continuous Closure in Abdominal Wound Repair: A Meta-analysis of 23 Trials. July 2008 ASIAN JOURNAL OF SURGERY VOL 31,NO 3, P104-114. 6. Chan V, Chan P-K, Chiu K-Y, Yan C-H, Ng F-Y. Does Barbed Suture Lower Cost and Improve Outcome in Total Knee Arthroplasty? A Randomized Controlled Trial. The Journal of Arthroplasty. 2017; 32; 1474-1477.

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Referencias:

1. Sundaram K, Warren J, Klika A, Piuze N, Mont M, Krebs V. Barbed sutures reduce arthroscopy closure duration compared to interrupted conventional sutures for total knee arthroplasty: a randomized control trial. Musculoskelet Surg. 2020;1-7.
2. Zayed M, Fouda U, Elsetohy K, Zayed S, Hashem A, Youssef M. Barbed sutures versus conventional sutures for uterine closure at cesarean section; a randomized controlled trial. The journal of maternal-fetal & neonatal medicine. 2017;1-8.
3. Ethicon, AST-2012-0331. Tissue gapping under tension of porcine cadaveric skin incisions closed with Stratafix Spiral in comparison to Monocryl in both interrupted and continuous stitching patterns. October 2012. Data on File.
4. Ethicon, AST-2013-0056. Performance Testing of STRATAFIX Symmetric PDS™ Size 2-0 suture device for Tissue Holding Strength with Multiple Incision Defects to Measure Gapping. April 2013. Data on File.
5. Ethicon, 100326296: Time zero tissue holding - Competitive claims comparisons for STRATAFIX Knotless Tissue Control Devices vs various products. 2015. Data on File

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
Referencias: 1. Ethicon, 100326296 Time Zero Tissue Holding - Competitive Claims Comparisons for STRATAFIX™ Knotless Tissue Control Devices vs Various Products. May 2015. Data on File. (060117-220713) 2. Ethicon, AST-2011-0210. Study to evaluate the tissue holding performance at time zero of DOLFIN PDS™ PLUS barbed suture sizes 1 and 2-0 vs dyed PDS™ II Plus suture sizes 1 and 2-0 in a continuous stitch pattern—Project DOLFIN 11822. July 2011. Data on File. (060117-220713) 3. Ethicon, AST-2013-0603. Performance Testing of STRATAFIX™ SYMMETRIC PDS™ PLUS Size 0 & 1 Devices – Initiation Strength in Porcine Tissue. April 2014. Data on File. (060117-220713)

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
Referencias:

1. Rodeheaver G, Píneros-Fernandez A, Salopek LS, et al. Barbed sutures for wound closure: in vivo wound security, tissue compatibility and cosmesis measurements. In: Transactions from the 30th Annual Meeting of the Society for Biomaterials; Mount Laurel, NK; p232.
2. Eickmann T, Quane E. Total knee arthroplasty closure with barbed sutures. J Knee Surg. 2010;23(3):163-167.
3. Ethicon, 100025465 rev;4, STRATAFIX™ SYMETRIC PDSTM Plus Knotless Tissue Control Device, Antibacterial. Instructions for Use. Dec 2019. Data on File.
4. Ethicon, 100376782 rev;6, STRATAFIX™ SPIRAL MONOCRYL™ Plus Knotless Tissue Control Device, Antibacterial. Instructions for Use. Nov. 2020. Data on File.
5. Ethicon, 100379555 rev;6 STRATAFIX™ SPIRAL PDSTM Plus Knotless Tissue Control Device, Antibacterial. Instructions for Use. Nov 2020. Data on File.
6. Ethicon, LAB-100028661v7- STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device Instructions for Use. Nov 2020. Data on File.
7. Ethicon, LAB-100088239v3 - STRATAFIX™ Spiral PDO Knotless Tissue Control Device Instructions for Use. Feb 2015. Data on File.
8. Ethicon, LAB100028658v3 STRATAFIX™ Spiral PGA-PCL Knotless Tissue Control Device. Instructions for Use. Sept 2012. Data on File.
9. Ethicon, 100326296: Time zero tissue holding - Competitive claims comparisons for STRATAFIX Knotless Tissue Control Devices vs various products. 2015. Data on File.
10. Sundaram K, Warren J, Klika A, Piuze N, Mont M, Krebs V. Barbed sutures reduce arthrotomy closure duration compared to interrupted conventional sutures for total knee arthroplasty: a randomized control trial. Musculoskelet Surg. 2020;1-7.
11. Zayed M, Fouda U, Elsetchy K, Zayed S, Hashem A, Youssef M. Barbed sutures versus conventional sutures for uterine closure at cesarean section; a randomized controlled trial. The journal of maternal-fetal & neonatal medicine. 2017;1-8. (060070-221220).
12. Ethicon, AST-2012-0331. Tissue gapping under tension of porcine cadaveric skin incisions closed with Stratafix Spiral in comparison to Monocryl in both interrupted and continuous stitching patterns. Aug 2012. Data on File. (060070-221220).
13. Ethicon, AST-2013-0056. Performance Testing of STRATAFIX Symmetric PDS Size2-0 suture devices for Tissue Holding Strength with Multiple Incision Defects to Measure Gapping. April 2013. Data on File. (060070-221220)

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

Referencia 1. Ethicon, 100326296 Time Zero Tissue Holding - Competitive Claims Comparisons for STRATAFIX™ Knotless Tissue Control Devices vs Various Products. May 2015. Data on File. (119236-230214)

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
Referencia:

1. Ethicon, 100326296 Time Zero Tissue Holding - Competitive Claims Comparisons for STRATAFIX™ Knotless Tissue Control Devices vs Various Products. May 2015. Data on File.

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



Referencia: 1. Ethicon LAB-100025466(4) -STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device Instructions for Use, Dec 2019, Data on File.

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* Triclosan has in vitro activity that inhibits bacterial colonization of the suture. For illustration purposes onl

Referencias: 1. Ming X, Rothenburger S, Nichols MM. In vivo and in vitro antibacterial efficacy of PDS Plus (polydioxanone with triclosan) suture. Surg Infect (Larchmt). 2008;9(4):451-457. (EM_ETH_WOUN_102315)

2. Ming X, Rothenburger S, Yang D. In vitro antibacterial efficacy of Monocryl Plus Antibacterial Suture (polygelcaprone 25 with triclosan). Surg Infect (Larchmt). 2007;8(2):201-207 (EM_ETH_WOUN_102315) 3. Bhende S, Burkley D, Nawrocki J. In vivo and in vitro antibacterial efficacy of absorbable barbed polydioxanone monofilament tissue control device with triclosan. Surg Infect 2018; Volume 19 (4):430-437. (EM_ETH_WOUN_102316)

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Global Health Authorities

World Health organization 2016:⁴
"The panel suggests the use of triclosan-coated sutures for the purpose of reducing the risk of SSI, independent of the type of surgery."
(Conditional recommendation, moderate quality of evidence)

Center for disease control and prevention, 2017:⁵
"Consider the use of triclosan-coated sutures for the prevention of SSI."

American College of surgeons/ Surgical infection society, 2016 SSI-guidelines:⁶
"Triclosan antibacterial sutures use is recommended for wound closure in clean and clean contaminated abdominal cases."

Local Health Authorities (UK and Germany)

NICE Guidelines Update 2019:⁷
"When using sutures, consider using antimicrobial triclosan-coated sutures, especially for pediatric surgery, to reduce the risk of surgical site infection (2019)."

Commission for Hospital Hygiene and Infection Prevention 2018:⁸
"Antiseptically coated sutures only reduce the threat of infection where there are high baseline SSI rates, contamination class III and IV surgeries, and patients with multiple morbidities (Cat.II)."

European Network

European Network for HTA 2017:⁹
"A statistically significant benefit of triclosan-coated sutures in reducing the risk of total incisional SSI was demonstrated in our systematic review /meta-analysis, based on moderate quality RCTs Data."

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
Referencias. 1. Ming X, Rothenburger S, Nichols M. In vivo and in vitro antibacterial efficacy of PDS Plus (polidioxanone with triclosan) suture. Surg Infect (Larchmt). 2008;9(4):451-457. (119401-230215) 2. Ming X, Rothenburger S, Yang D. In vitro antibacterial efficacy of Monocryl Plus Antibacterial Suture (poliglecaprone 25 with triclosan). Surg Infect (Larchmt). 2007;8(2):201-207. (119401-230215) 3. Barbolt T. Chemistry and safety of triclosan, and its use as an antimicrobial on Coated VICRYL* Plus Antibacterial Suture (coated polyglactin 910 suture with triclosan). Surg Infect (Larchmt). 2002;3(suppl1):S45-S53. (119401-230215) 4. World Health Organization. Global guidelines for the prevention of surgical site infection. <https://apps.who.int/iris/bitstream/handle/10665/250680/9789241549882-eng.pdf?sequence=8>. Published November 2016. Accessed July 13, 2022. (156562-230215) 5. Berrios-Torres S, Umscheid C, Bratzler D, et al. Centers for disease control and prevention guideline for the prevention of surgical site infection, 2017. JAMA Surg. 2017;152(8):784-791. (156562-230215) 6. Ban K, Minei J, Laronga C, et al. American college of surgeons and surgical infection society: surgical site infection guidelines, 2016 update. J Am Coll Surg. 2017; 224(1):59-74. (156562-230215) 7. NICE Guideline Updates Team (UK). Surgical site infection: prevention and treatment. NICE website. <https://www.nice.org.uk/guidance/ng125/chapter/Recommendations#closuremethods> Accessed July 13, 2022. (156562-230215) 8. Prevention of postoperative wound infections. Recommendation of the Committee for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute. Bundesgesundheitsbl. 2018;61(4):448-473. (156562-230215) 9. EUnetHTA. Rapid assessment of JA3 WP4 on "Antibacterial-Coated Sutures Versus Non-Antibacterial-Coated Sutures For The Prevention Of Abdominal, Superficial And Deep Incisional, Surgical Site Infection (SSI) March 2017, V1.4. (156562-230215)

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Referencias: 1. © NICE 2021. MEDICAL TECHNOLOGY GUIDANCE: PLUS SUTURES FOR PREVENTING SURGICAL SITE INFECTION. Available from: nice.org.uk/guidance/MTG59 Accessed on: 28 June 2021. All rights reserved. Subject to Notice of rights. NICE guidance is prepared for the National Health Service in England. All NICE guidance is subject to regular review and may be up dated or withdrawn. NICE accepts no responsibility for the use of its content in this product/publication. (187177-230216)

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




Claramente identificables

STRATAFIX™ Spiral
Unidireccional

STRATAFIX™ Spiral Device
Bidirectional

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Referencias: 1. Ethicon, 389974R06 STRATAFIX™ Spiral PDS™ Plus Unidirectional Knotless Tissue Control Device, Instructions for Use, November 2020. Data on File. 2. Ethicon, 389972R06 STRATAFIX™ Spiral MONOCRYL™ Plus Unidirectional Knotless Tissue Control Device, Instructions for Use, November 2020. Data on File. 3. Ethicon, 489974R01 STRATAFIX™ Spiral PDS™ Plus Bidirectional Knotless Tissue Control Device, Instructions for Use, December 2021. Data on File. 4. Ethicon, 389947R01 STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Knotless Tissue Control Device, Instructions for Use, December 2021. Data on File. 5. Ethicon, 395948R01 STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device, Instructions for Use, December 2021. Data on File. 6. Ethicon, 03-5437R5 STRATAFIX™ Spiral PDO Knotless Tissue control Device, Instructions for Use, February 2021. Data on File. 7. Ethicon, 03-5438R5 STRATAFIX™ Spiral PGA-PCL Knotless Tissue control Device, Instructions for Use, February 2021. Data on File. 8. Ethicon, 03-5439R5 STRATAFIX™ Spiral Polypropylene Knotless Tissue control Device, Instructions for Use, February 2021. Data on File. 9. Ethicon, 03-5445R4 STRATAFIX™ Spiral Polypropylene Unidirectional Knotless Tissue control Device, Instructions for Use, February 2021. Data on File.

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Johnson & Johnson, S.A.
© Johnson & Johnson S.A. Todos los derechos reservados